



AG connectors are made from a polypropylene polymer graded for biopharmaceutical industries.

Light, easy to use and designed to guarantee high level of security for your Single-Use assemblies.

TYPICAL APPLICATIONS

- Fluid and media transport in biopharmaceutical processes and biotechnology
- Single-use assemblies

FEATURES

- Manufactured in ISO Class 7 cleanroom
- Sterilization Methods: Autoclave steam, Gamma radiation, 25-45 kGy
- Regulatory data package
 - FDA 21 CFR 177.1520
 - Polypropylene Resin used for manufacturing Tri-clover gasket is successfully tested for USP Class VI (Toxicity test) from independent 3rd party lab.
 - USP 87 certified
 - ISO 10993
 - TSE/BSE free

IMPORTANT: It is the user's responsibility to ensure the suitability and safety of Aseptic Group products for all intended uses. Products must be used in accordance with applicable regulatory requirements defined in datasheet. Aseptic Group assumes no obligations or liability for any advice furnished by it, or for results obtained with respect to these products. The data supplied is given and accepted at the buyer's risk.

WARRANTY: Aseptic Group products are warranted to be free from defects in materials and workmanship. Any products found to carry such defects will be replaced without charge. This warranty is limited only to the replacement of Aseptic Group products. Aseptic Group specifically disclaims any other liability.



BENEFITS

- Smooth transition design
- Eliminates entrapment and leak points that occur with traditional fitting systems
- No excessive voids at the interface or in the fluid path
- Compatible with flexible tubing
- PP is gamma stable

09/2022



PRODUCT SPECIFICATIONS

AG-TC075	
Part Number	For tubing dimensions inches (mm)
AG-TC075-18-PP	1/8 (3,2)
AG-TC075-14-PP	1/4 (6,4)
AG-TC075-38-PP	3/8 (9,5)
AG-TC075-12-PP	1/2 (12,7)
AG-TC075-34-PP	3/4 (19,1)

AG-TC150	
Part Number	For tubing dimensions inches (mm)
AG-TC150-14-PP	1/4 (6,4)
AG-TC150-38-PP	3/8 (9,5)
AG-TC150-12-PP	1/2 (12,7)
AG-TC150-34-PP	3/4 (19,1)
AG-TC150-1-PP	1 (25,4)

AG-WYE	
Part Number	For tubing dimensions inches (mm)
AG-WYE-14-PP	1/4 (6,4)
AG-WYE-38-PP	3/8 (9,5)
AG-WYE-12-PP	1/2 (12,7)
AG-WYE-34-PP	3/4 (19,1)

AG-TEE	
Part Number	For tubing dimensions inches (mm)
AG-TEE-14-PP	1/4 (6,4)
AG-TEE-38-PP	3/8 (9,5)
AG-TEE-12-PP	1/2 (12,7)
AG-TEE-34-PP	3/4 (19,1)



IMPORTANT: It is the user's responsibility to ensure the suitability and safety of Aseptic Group products for all intended uses. Products must be used in accordance with applicable regulatory requirements defined in datasheet. Aseptic Group assumes no obligations or liability for any advice furnished by it, or for results obtained with respect to these products. The data supplied is given and accepted at the buyer's risk.

WARRANTY: Aseptic Group products are warranted to be free from defects in materials and workmanship. Any products found to carry such defects will be replaced without charge. This warranty is limited only to the replacement of Aseptic Group products. Aseptic Group specifically disclaims any other liability.

09/2022

TUBINGS

QUICK
COUPLERS

STERILE CONNECTIONS AND DISCONNECTIONS /
SEALER SYSTEMS AND TUBING WELDERS

GASKETS

ADAPTERS
AND FITTINGS

CONTAINERS

CLOSURES

SENSORS AND
MONITORS

PUMPS




PRODUCT SPECIFICATIONS

AG-X	
Part Number	For tubing dimensions inches (mm)
AG-X-14-PP	1/4 (6,4)
AG-X-38-PP	3/8 (9,5)
AG-X-12-PP	1/2 (12,7)
AG-X-34-PP	3/4 (19,1)

AG-C	
Part Number	For tubing dimensions inches (mm)
AG-SC-12-PP	1/2 (12,7)
AG-SC-14-PP	1/4 (6,4)
AG-SC-38-PP	3/8 (9,5)
AG-SC-34-PP	3/4 (19,1)

AG-RC	
Part Number	For tubing dimensions inches (mm)
AG-RC-14-18-PP	1/4 (6,3) to 1/8 (3,2)
AG-RC-12-18-PP	1/2 (12,7) to 1/8 (3,2)
AG-RC-38-14-PP	3/8 (9,5) to 1/4 (6,4)
AG-RC-34-14-PP	3/4 (19,1) to 1/4 (6,4)
AG-RC-12-38-PP	1/2 (12,7) to 3/8 (9,5)
AG-RC-34-38-PP	3/4 (19,1) to 3/8 (9,5)

AG-PLUG	
Part Number	For tubing dimensions inches (mm)
AG-PLUG-18-PP	1/8 (3,2)
AG-PLUG-14-PP	1/4 (6,4)
AG-PLUG-38-PP	3/8 (9,5)
AG-PLUG-12-PP	1/2 (12,7)
AG-PLUG-34-PP	3/4 (19)



IMPORTANT: It is the user's responsibility to ensure the suitability and safety of Aseptic Group products for all intended uses. Products must be used in accordance with applicable regulatory requirements defined in datasheet. Aseptic Group assumes no obligations or liability for any advice furnished by it, or for results obtained with respect to these products. The data supplied is given and accepted at the buyer's risk.

WARRANTY: Aseptic Group products are warranted to be free from defects in materials and workmanship. Any products found to carry such defects will be replaced without charge. This warranty is limited only to the replacement of Aseptic Group products. Aseptic Group specifically disclaims any other liability.

09/2022

TUBINGS

QUICK
COUPLERS

STERILE CONNECTIONS AND DISCONNECTIONS /
SEALER SYSTEMS AND TUBING WELDERS

GASKETS

ADAPTERS
AND FITTINGS

CONTAINERS

CLOSURES

SENSORS AND
MONITORS

PUMPS